

Rapid

F217/A/B/C

Oil Extraction Apparatus

To determine the oil or finish content of wool or synthetic samples by solvent percolation and evaporation. Suitable for fiber, yarn or fabric. Available with 1, 2 or 3 extraction units. Complies with Edana method ERT 155(needs S217014 large extraction tube fully complaint). Full title of new method is: DETERMINATION OF ACETONE EXTRAC - TABLE FINISH ON NONWOVENS.

Specifications	F217/A/B/C
Weight	4.4 kg
Dimensions	170 x 170 x 380 mm (WxDxH)
Standards	Edana method ERT 155



F217/A/B/C

- To determine the oil or finish content of wool or synthetic samples by solvent percolation and evaporation.
- New method: Determination Of Acetone Extrac - Table Finish On Nonwovens.